EN  
3

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













Rekommenderad bruksanvisning

# The Vital-Port® Vascular Access System

## INTRODUCTION

The Vital-Port® Vascular Access System is a completely implantable subcutaneous system for patients who require repeated intermittent vascular access for simultaneous or separate infusion and/or withdrawal of fluids. The Vital-Port® System consists of either one or two independent injection port chambers, each with a self-sealing septum accessed by percutaneous needle puncture using a non-coring needle, and a catheter.

## Product Description

The Vital-Port® Vascular Access System is supplied sterile and non-pyrogenic. The titanium or polysulfone port body has suture points that permit anchoring. These are silicone filled in titanium models. A silicone or polyurethane catheter is supplied with a catheter locking sleeve. The catheter may be supplied pre-attached or detached. A non-coring stainless steel needle is included with all models. A vascular access introducer set, tunneler and a non-coring infusion set may also be included with the Vital-Port® system. All Vital-Port® systems are for single use only **Figure 1** A=Height B=Length C=Width D=Diameter E=Suture Points F=Septum G=Catheter H=Catheter Lock.

## Indications for Use

The Vital-Port® Vascular Access System is indicated for patient therapy requiring repeated vascular access for simultaneous or separate injection or infusion therapy and/or blood sampling.

## Contraindications

The Vital-Port® System should not be implanted in the presence of known or suspected infections, bacteremia, septicemia, or peritonitis, or in patients who have demonstrated prior intolerance to the materials of construction.

## Potential Complications

Use of the Vital-Port® System involves potential risks normally associated with the insertion or use of any implanted device or indwelling catheter, including but not limited to: air embolism; bacteremia, infection, sepsis; brachial plexus injury; cardiac arrhythmia, puncture, tamponade; catheter disconnection, fracture, fragmentation, malposition, pinch-off, rupture, shearing; embolization, extravasation, formation of calcium-like deposits; fibrin sheath formation; hemothorax; hemoatoma; port/catheter erosion, migration, occlusion, rejection, rotation; pneumothorax; thoracic duct injury; thrombosis; thromboembolism; vessel erosion, laceration, perforation, trauma.

## WARNINGS

• Introduction of the catheter into the subclavian vein **Figure 2A** using standard percutaneous techniques may subject the catheter to periodic compression forces within the narrow costoclavicular space between the clavicle **Figure 2B** and first rib **Figure 2C**. Reported complications from repeated subclavian compression include catheter pinch-off syndrome, catheter fracture, and catheter shear followed by embolization of the distal portion. The likelihood of catheter occlusion and damage can be greatly reduced by choosing an access site lateral to the clavicle (preferably at or lateral to the midclavicular line) **Figure 2D** and prior to the vein entering the thorax at the costoclavicular space.


RECOMMENDED: Central Venous Catheter introduced at or lateral to midclavicular line **Figure 2E**.

AVOIDED: Central Venous Catheter introduced medial to costoclavicular space **Figure 2F**.

- Prevent air embolism by applying finger pressure over the opening of the introducer sheath.
- Failure to adequately anchor the port to the fascia increases the risk of catheter fracture and/or disconnection which could result in catheter migration.

## Peripheral Implanted Port

- To avoid catheter damage to a vascular access system implanted peripherally in the upper arm, position the catheter so that it enters the venous system in the lower third of the upper arm and below the vein's passage through the deep or brachial fascia. Following insertion, confirm that the catheter follows a smooth course into the vein and that no abrupt directional changes are present.

- Abrupt directional changes may result in catheter fracture and subsequent possible embolization. A radiographic confirmation of the catheter insertion should be made to ensure that the catheter is not being pinched.

## PRECAUTIONS

### General

- Consult the “Suggested Instructions for Use” before implanting or accessing the product.
- Do not use this product if the expiration date has passed. Do not use this product if the package appears to be opened or damaged.
- Use aseptic technique when implanting or accessing the device.
- Use talc-free gloves for handling the device.
- To avoid the potential problems associated with catheter damage incurred due to unnecessary or improperly maintained long-term device placement, consider removing the entire port system following cessation of therapy. A post-operative x-ray is suggested for verification of system removal.

### Implant Considerations

- Withdrawal or manipulation of distal spring coil portion of wire guide through needle tip may result in breakage.
- Avoid positioning the port septum directly below the incision line since the septum will be punctured repeatedly.
- Avoid contacting the catheter with sharp objects which may lead to catheter shearing.
- Avoid excessive clamp force which could damage the catheter. Clamp the catheter only on a section that will be cut off prior to implant. If the catheter was previously cut to length, keep the distal end elevated rather than clamping.
- Avoid trimming the catheter at an angle: an angled catheter tip can be sealed off against the side of a vessel.
- Improper catheter connection may result in catheter damage, leakage, or possible disconnection.
- Do not practice making a port-catheter connection with a catheter intended for implantation as damage to the catheter may result. Prior practice is recommended to ensure dexterity in connecting the catheter to the port.
- Verify that the catheter lock and reinforcing sleeve are retrieved along with the port body and total length of the catheter upon removal of the Vital-Port® System. A post operative x-ray is suggested for verification.

### Accessing the Port

- Follow all manufacturer guidelines for infusates.
- Do not exert excessive force during needle insertion. The needle point may become damaged which could tear the septum when the needle is removed.
- Do not exceed 40 psi pressure when administering fluid through the system. Use of a 10 ml or larger syringe is recommended for all injection and flushing procedures. Pressures in excess of 40 psi can be generated with most syringes: the smaller the syringe the higher the pressure that can be generated. For VITAL-PORT® vascular access system models that are labeled as a POWER INJECTABLE PORT, refer to the *Vital-Port Power Injection System Suggested Instructions for Use* that is provided with the product.
- Inspect and replace the infusion set and sterile dressing according to CDC guidelines, or at intervals according to hospital policy or current practice.

## PATIENT COUNSELING INFORMATION

### Patient Resources

Each Vital-Port® System is supplied with a patient information card, containing the port model and lot number as well as manufacturer contact information. Please complete the requested information on the card and instruct the patient to carry the card with them.

A “Patient Information” booklet is also supplied with each port or available separately.

## IMPLANTATION PROCEDURE FOR CENTRAL VENOUS VITAL-PORT® SYSTEMS

### System Preparation

Prior to placement, all air must be flushed from the chamber(s) and catheter lumen(s) of the Vital-Port® system using a non-coring needle, a 10ml or larger syringe and heparinized saline (100 IU/ml).

**CAUTION:** Use aseptic techniques and extreme care when handling the catheter. Avoid both contacting the catheter with sharp objects, which may lead to catheter shearing, and excessive clamp force, which could damage the catheter. Clamp the catheter only on a section that will be cut off prior to implant. If the catheter was previously cut to length, keep the distal end elevated rather than clamping. When trimming catheter to size, perform perpendicular cut in one motion.

### Port Placement Considerations

The injection port reservoir should be placed where it will be supported by an underlying bony structure. The location should be convenient and comfortable for the patient. Patient mobility should also be considered.

The incision should not be significantly larger than the length of the port base. Avoid positioning the port septum directly below the incision line since the septum will be punctured repeatedly.

The port should not be placed too deeply to be easily located, or too superficially, causing erosion due to pressure necrosis. A depth between 0.5cm and 2cm under the skin surface is recommended. The port should be easily located by palpation for access.

### Catheter Placement Considerations

The catheter tip should be placed in an area of high blood flow (usually the SVC) **Figure 3**.

**A=Axillary Vein B=Subclavian Vein C=Innominate Vein D=Superior Vena Cava**  
**E=Catheter Tip Positioning zone.**

### WARNING

• Introduction of the catheter into the subclavian vein **Figure 2A** using standard percutaneous techniques may subject the catheter to periodic compression forces within the narrow costoclavicular space between the clavicle **Figure 2B** and first rib **Figure 2C**. Reported complications from repeated subclavian compression include catheter pinch-off syndrome, catheter fracture, and catheter shear followed by embolization of the distal portion. The likelihood of catheter occlusion and damage can be greatly reduced by choosing an access site lateral to the clavicle (preferably at or lateral to the midclavicular line) **Figure 2D** and prior to the vein entering the thorax at the costoclavicular space.

**RECOMMENDED:** Central Venous Catheter introduced at or lateral to mid-clavicular line **Figure 2E**.

**AVOIDED:** Central Venous Catheter introduced medial to costoclavicular space **Figure 2F**.

For additional information and references on this subject, a booklet entitled Understanding and Avoiding Subclavian Compression of Central Venous Catheters is available from Cook.

### Accessing The Vein

To allow for body movement, the catheter should not be implanted with tension. Also avoid excessive catheter slack in the port pocket to prevent possible kinking.

### Percutaneous Procedure

- Use an appropriately sized introduction system.
- Attach introducer needle to syringe.
- Introduce needle into desired vessel, gently aspirating during introduction.
- Remove syringe from needle being sure to cover needle opening to prevent air embolism.

**WARNING: Prevent air embolism by applying finger pressure over the opening of the introducer sheath.**

- Insert wire guide into needle using wire straightener.
- Advance the wire guide to the appropriate position. Confirm using radiographic technology.
- Remove needle and straightener.
- Advance introducer dilator/sheath over wire guide.
- Remove dilator.
- Insert catheter into sheath. Position distal end of catheter at desired location using applicable technologies.

- k. Peel away sheath while withdrawing from vessel.
- l. Verify correct catheter tip placement using fluoroscopy or other appropriate technology.

#### **Cut-down procedure**

- a. Make a small incision to expose the entry vessel of choice.
- b. Isolate and stabilize the vessel. Perform vessel incision.
- c. Introduce catheter and advance tip to desired location.
- d. Secure catheter at venotomy site.
- e. Verify correct catheter tip placement using fluoroscopy or other appropriate technology.

#### **Pre-attached Catheter Models**

1. Prior to catheter placement, select the site for the port pocket.
2. Measure the appropriate catheter length. Provide enough slack in the catheter to allow for body movement.
3. Trim the excess by cutting the distal end squarely.
4. Verify that the port chamber(s) and catheter are filled with heparinized saline (100 IU/ml).
5. Following a standard percutaneous procedure or surgical cutdown procedure, gain access to the vascular system. Use a standard tunneling device when required.
6. Verify catheter positioning radiographically.

#### **Detached Catheter Models**

1. Verify that the catheter is filled with heparinized saline (100 IU/ml).
2. Gain access to the vascular system using standard percutaneous techniques or surgical cutdown.
3. Position the distal end of the catheter in the vessel at the desired location. Advance the catheter from the selected access vessel back to the port pocket. Use a standard tunneling device when required.
4. After placement, confirm correct positioning of the catheter using a radiographic techniques.
5. To allow for body movement, the catheter should not be implanted with tension, but avoid excessive catheter slack in the port pocket to prevent possible kinking.

**CAUTION: Trim the proximal end squarely to aid proper sealing of the locking mechanism.**

6. Verify that the port chamber(s) and outlet tube(s) are filled with heparinized saline (100 IU/ml), then attach the catheter to the port outlet tube(s) as indicated.

#### **Peripheral Models**

The Vital-Port® Mini model 5116W is supplied with a detached catheter having a Slip-Coat™ Wire Guide Obturator pre-inserted for facilitating peripheral placement.

#### **PRECAUTIONS**

The Vital-Port® Mini peripheral system should be used only after careful consideration and with extreme caution for patients who have had extensive trauma to the veins of their arms or have anatomical irregularities which may hamper proper placement of the system.

1. Connect a syringe with heparinized saline (100 IU/ml) to the stopcock and flush the catheter. Close the stopcock by rotating the handle 90° **Figure 4A**.
2. The catheter should be placed via a percutaneous puncture technique or a cutdown procedure.
3. Position the distal end of the catheter in the vessel at the desired location.
4. With the catheter in the desired position, remove the obturator from the catheter by removing the stopcock **Figure 5A** and gently withdrawing the obturator **Figure 5B**.
5. Next, flush the catheter with heparinized saline (100 IU/ml) via the Luer connector on the end of the catheter.
6. Clamp the catheter closed at the proximal end to be connected to the portal, clamping the catheter only on the section that will be cut off prior to implant. Cut the Luer connector off of the catheter.

#### **SYSTEM ASSEMBLY**

##### **Single Chamber Models**

1. Position the catheter locking sleeve over the proximal end of the catheter **Figures 6A, 7A** and push the catheter over the port outlet **Figures 6B and 7B**. To facilitate advancement over the ring, rotate the catheter while pushing forward.
- A. On Standard size models, the end of the catheter should be advanced halfway **Figure 6D** between the ring **Figure 6C** and the shoulder of the connector tube **Figure 6E**.

8. With Petite and Mini models, the end of the catheter should be advanced halfway **Figure 7D** between the ring **Figure 7C** and port body **Figure 7E**.
9. Slide the catheter locking sleeve over the port outlet tube and catheter until it is engaged in the recess of the port body.
10. When properly positioned, the end of the locking sleeve should not be visible. The end of the catheter should be visible through the clear portion of the catheter locking sleeve **Figure 8**.

**CAUTION: Improper catheter connection may result in catheter damage, leakage, or possible disconnection. Prior practice is recommended to ensure dexterity in connecting the catheter to the port. Do not practice with a catheter intended for implantation as damage to the catheter may result.**

#### **Dual Chamber Models**

1. Position the catheter locking sleeve over the proximal end of the catheter and push the catheter over the port outlet tubes, **Figure 9A** aligning the catheter lumens with the port outlet tubes.
2. The end of the catheter should be positioned against the port body, completely covering the port outlet tubes. This can be visualized through the window in the bottom of the port body **Figure 9B**.
3. Slide the catheter locking sleeve over the port outlet tubes and catheter until engaging the threads in the port body. Turn the catheter locking sleeve clockwise until rotation is stopped by the locator pin. When properly positioned, the wings of the catheter locking sleeve will be parallel with the base of the port.
4. The end of the catheter should now appear as a small "O" ring or donut within the port body window **Figure 9C**.

**CAUTION: Improper catheter connection may result in catheter damage, leakage, or possible disconnection. Prior practice is recommended to ensure dexterity in connection the catheter to the port. Do not practice with a catheter intended for implantation as damage to the catheter may result.**

#### **Flow Verification**

Before closing the port pocket, unimpeded fluid flow should be verified. To confirm two-way fluid flow, use aspiration followed by a flush with 20 ml saline to cleanse the catheter lumen(s) and port chamber(s) of blood. Flush each chamber and corresponding catheter lumen with 5 ml of heparinized saline (100 IU/ml) (3 ml for Mini systems). For dual models, check both chambers and corresponding catheter lumens. Use a non-coring needle and a 10ml or larger syringe with saline to confirm that the flow is not obstructed and that there are no leaks. After flow verification, establish a heparin lock with 5ml of heparinized saline (100 IU/ml) (3 ml for Mini systems). Continue to inject fluid and depress the port with the fingers while withdrawing the needle to avoid reflux.

#### **Port Sutures and Site Closure**

Anchor the port to the fascia with nonabsorbable sutures. Three sutures, with at least one at each end of the port, are recommended. Close the site, making sure that the port does not lie directly beneath the incision. Access the port, to ensure a proper heparinized saline lock, before sending the patient to the recovery area.

**WARNING: Failure to adequately anchor the port to the fascia increases the risk of catheter fracture and/or disconnection which could result in catheter migration.**

#### **PROCEDURES FOR USING THE VITAL-PORT® SYSTEM**

##### **Accessing the Port Chamber**

- Always use aseptic technique before accessing the system.
- Locate the port by palpation.
- Observe for wound hematoma or accumulation of serous fluid around the implant site.

**Syringes:** Use 10 ml syringes or larger for all flushing or injection procedures.

**Needles:** Use only non-coring needles, 19 to 22 gauge. Non-coring needles preserve the life of the self-sealing septum. Needles are for single use only.

**Needle Insertion:** Needles should be inserted perpendicular to the implanted port after locating the port septum. Puncture the skin immediately over the septum and advance the needle point slowly through the septum until it makes contact with the bottom of the port chamber **Figure 10**. Incomplete needle insertion may result in needle obstruction by the septum. Never leave the system open to air. Do not rock or tilt the needle. Stabilize the needle in the port septum and cover with sterile dressing.

**Saline Flushes:** Prior to drug administration, flush the system with saline solution. If more than one drug is administered during the treatment, flush the system with saline solution between drugs. After completing patient treatment, always flush the system with saline to cleanse the catheter and port chamber.

**Heparin Lock:** To keep the Vital-Port® system patent, establish a heparin lock by flushing each port chamber and catheter lumen with 5 ml heparinized saline (100 IU/ml) (3 ml for Mini systems) after each patient treatment and at least once every 4 weeks.

### **Injection/Infusion**

Before accessing the port chamber, connect the syringe or IV set to a non-coring needle, and fill the infusion set with saline to remove air.

Access the port as previously described and begin infusion or injection. Examine the puncture site closely for any signs of subcutaneous irritation.

Flush the system as previously described.

### **Blood Sampling**

Flush the port with saline to confirm system patency. Withdraw 3-4 ml blood sample and discard. Withdraw the required blood volume. Immediately follow the sampling with a 20 ml saline flush to cleanse the port chamber and catheter. If no further treatments are required, establish a heparin lock as previously described.

### **DISCONTINUING USE OF THE VITAL-PORT® SYSTEM**

In vitro testing has shown sealing integrity after 1200 and 2500 punctures for 20 and 22 ga, respectively (400 and 1000 punctures for the Mini). Septum lifetime may vary, depending on clinical use.

It is suggested that the clinician consider removing the Vital-Port® System upon completion of treatment.

Upon removal of the Vital-Port® System, verify that the catheter lock and reinforcing sleeve are retrieved along with the port body and total length of the catheter. A post operative x-ray is suggested for verification.

Ports and catheters that are removed because the clinician believes there has been a system malfunction should be returned to Cook for evaluation. Please contact customer service for a returned goods authorization number and instructions.

### **HOW SUPPLIED**

This device is supplied sterile and intended for one-time use.

If device packaging is damaged, do not use and discard the device.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

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#### EN MR-Conditional

Testing has demonstrated that the Vital-Port™ Vascular Access Systems are MR-conditional. A patient with the implant can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less • Spatial gradient magnetic field of 720 Gauss/cm or less
- Maximum MR-system reported whole body-average specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning

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Vital-Port® is a registered trademark of COOK Vascular™ Incorporated. COOK® is a registered trademark of Cook Incorporated. Slip-Coat™ is a trademark of Sterilization Technical Services.